

Cardio Track 3C

Electrocardiograph



Cardio Track 3C

Electrocardiograph



Automatic ECG waveform and parameter measurement

Complete digital design, wave digital filter. Automatic adjustment of baseline, amplification, and automatic lead switch

Single Color TFT Display (Resolution: 320×240). Thermal Recorder (80mm)

Simultaneously 12 Leads Acquisition. Optional Selections for Rhythm Lead

High Resolution Thermal Printer, Synchronous Record Details of ECG Waveforms, name of leads, paper feed speeds, gain, patients' information, and analysis reports, etc.

Multi-language operation interface, Various Printout Format

Built-in Li-ion Battery, Support Uninterrupted Working Up to 2 Hours

150 Digital ECG Records

Optional Network Capability



Light weight and portable, AC/DC power supply



Simple interface display, easy operation mode.



Automatic measurement report



Excellent quality CE approved