



Fetal Monitor

Model : Betaview Minimax

Features

- Light weight, compact, table top or wall mount fetal monitor
- 5.7" LCD display screen and angle of the screen can be adjusted from 0° to 60°
- Basic Functions : FHR, TOCO, Manual fetal movement marking
- Highly sensitive broad band pulsed wave transducer with advance autocorrection & DSP techniques.
- Visual & audible alarm for upper & lower limit FHR, TOCO
- Built-in thermal recorder.
- Built-in RS-485 network interface and rechargeable battery backup for 4 hours.

Specification

Ultrasound (FHR)

ULTRASOUND FREQUENCY	1.0MHZ Pulsed wave transducer
FHR RANGE	65bpm - 210bpm
PRECISION	±0.25bpm
RECORD DIFF.	30bpm/cm

UA

RANGE	0 - 100units
PRECISION	±1unit

Event Marker

Manual Fetal Movement Event Marking

Recorder

PAPER SPEED 1,2,3 cms/min (optional)

Rechargeable Battery

BATTERY TYPE	Lithium-ion Rechargeable Battery
APPRO. CHARGING	12 Hours
BATTERY BACKUP	4 Hours

Optional : Upgradable to Maternal Monitoring
(ECG, NIBP, SpO2, TEMP).

Regd. Office

Vin BioTech Systems Limited

No. 2, II Cross Street, M.E.S. Street,
Tambaram East, Chennai - 600 059.
Ph : + 91 - 44 - 2239 1239, 6546 2737.
Email : sales@vinbiotech.com
Web Site : www.vinbiotech.com



ISO 9001 - 2000