

*Trusted
Your Partner in Healthcare*

Nebutech

(Air Compressing Nebulizer)



NEBUTECH COMPRESSOR NEBULIZER FEATURES:

1. Heavy duty ball bearing piston for use in hospital, clinics and home use.
2. Small size for convenient carrying
3. Reliable pressure and air flow
4. Continuous treatment > 5 hours (under room temperature)
5. Build-in carrying handle

MOTOR SPECIFICATIONS:

Motor	: Oil-free ball bearing piston. 6330 pump
Power: AC 230V or 120V	: 50 Hz/60Hz
Power consumption	: < 72W
Compressor Pressure range	: 30-50 psi; 0.21- 0.35 Mpa
Operation Pressure Range	: 12-25psi; 0.10-0.15 Mpa
Compressor Operation flow	: 12-16 LPM

PRODUCT PACKAGE:

Color Box	: Master Carton
Qty / ctn	: 1/5units
Net Weight	: 10.5 kgs
Gross Weight	: 12 kgS
Compressor Operation flow	: 12-16 LPM
Dimension	: 297 *180*108mm
Weight	: 1.8kgs



VIN BIOTECH SYSTEMS LIMITED

R.O.: No. 2, II Cross Street, M.E.S. Road, Tambaram (E), Chennai - 59.
W: No. 2, U.V. Swaminathan Iyer Street, M.E.S. Road, Tambaram (E),
Chennai - 600 059.
Tel: +91 - 44 - 2239 0795, 6456 2736
Email: info.vbsl@gmail.com, ceo.vinbiotech@gmail.com,
Website: www.vinbiotech.com, Skype: shani.p7

