

EDAN



SE-1200 series Twelve-channel ECG

User-friendly Design

- 8" high resolution color touch screen (SE-1200 Express)
- 5.7" LCD screen with blue backlight (SE-1200)
- Alphanumeric keyboard

Streamlined Workflow

- One-touch operation
- Support barcode scanner to simplify workflow

Accurate Clinical Performances

- Weak signal detection with the most advanced 24 bits A/D converter technology
- Anti-noise technology with comprehensive filters and the high CMRR over 115dB
- 120s real-time ECG waveforms freezing and review to assist your diagnosis
- Automatic measurements and interpretation tested by CSE/AHA/MIT database
- Extend printing once arrhythmia is detected to further rhythm analysis

Flexible Printing Solutions

- Built-in high resolution thermal recorder to print on normal thermal paper and A4 size fax paper
- External USB printer support

Data Management

- Large internal memories: 200 ECGs (SE-1200 Express)/100 ECGs (SE-1200), which can be enlarged by USB flash disk
- LAN/WIFI/RS232 connection to PC
- Smart ECG Viewer software for data management on PC (option)
- PDF, SCP and FDA-XML data export (option)



Touch Screen



Wireless AP



One-touch Operation

PC-based ECG Data Management Software



Care for Health